



## 離心濃縮組

### CENTRIFUGAL EVAPORATOR COMPLETE SET



#### 特點:

- 集離心濃縮機, 低溫冷凝機, 及真空幫浦於一組中, 不但經濟, 節省空間, 而且是絕配組合。
- 用於生化實驗 DNA, RNA, PEPTIDES之製備, 藥品分析及食品科學研究。
- 可達100%之樣品回收率。
- 整組合:
  1. CES-8000 離心濃縮機  
需另選配ROTOR
  2. CT-5020 (-50°C) 桌上型低溫冷凝機
    - 可選擇歧管座作為簡易型冷凍乾燥機
    - 溫度LED顯是為選擇配件
    - 可另選購-85°C機型。
  3. G-100D 油迴轉式真空幫浦
  4. 離心濃縮專用架

#### SPECIFICATIONS

MODEL	CES-8000
ROTOR SPEED	100~1,800 rpm with meter
TEMP. RANGE	R.T.+5°C~8°C
TEMP. STABILITY	±0.1°C
TEMP. SETTING & DISPLAY	LED display
DRIVE SYSTEM	Magnetically coupled
VACUUM GAUGE	0~760mmHg
CHAMBER MATERIAL	SUS304 plus teflon coating
CHAMBER DIMENSIONS	φ142 x 110mm H
ROTOR	Optional
TOTAL DIMENSIONS (WxDxH)	330 x 300 x 290mm
POWER SOURCE	110/220V, 50/60Hz
	CT-5020
TEMPERATURE	-50°C
CAPACITY	3.5ℓ (1ℓ/batch)
TRAP MATERIAL	SUS304 plus Teflon coating
SAFETY DEVICE	Over current prevention devices
POWER SOURCE	110/220V, 50/60Hz
DEFROST	Hot water
ACCESSORY	Ice-hand x 1
TOTAL DIMENSIONS (WxDxH)	340 x 450 x 490mm
VACUUM PUMP	6.7x10 <sup>-1</sup> Pa, 120ℓ/min.
STAND FOR COMPLETE SET	Standard

\*\* (本公司保有修改規格之權利, 不另行通知。)

We reserve the rights to change the specifications without notice.)